

<b>Product Name:</b>	<b>Serological Pipette</b>		
<b>Product Application:</b>	<b>NEST serological pipettes can be used for liquid transfer in biological experiments such as cell culture and bacterial culture.</b>		
<b>Product Materials:</b>	Polystyrene (PS).		
<b>product Specifications:</b>	Temperature range: stored at room temperature Shelf life: 3 years after date of production		
<b>Sterilization:</b>	<b>Yes.</b> Sterilized by E-beam, Sterility Assurance Level: SAL=10 <sup>-6</sup> . The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137		
<b>Pyrogens:</b>	Non-Pyrogenic		
<b>RNase/DNase:</b>	DNase/RNase free		
<b>BSE/TSE:</b>	These products are deemed animal free.		
<b>Performance Testing:</b>	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified                      Stress test: qualified Sealing test: qualified                                      Packaging inspection: qualified		
<b>Features</b>	<ol style="list-style-type: none"> <li>1. Made from high clarity polystyrene.</li> <li>2. Withstand temperature range: Approx. 0°C to 40°C</li> <li>2. Individual packaging.</li> <li>3. Edged color-coded band for identifying Pipettes of different volumes.</li> <li>4. Negative graduations for extra volume.</li> <li>5. Sterilized by E-beam.</li> <li>6. Non-pyrogenic, non-cytotoxic and non-hemolytic.</li> <li>7. DNase/RNase free.</li> <li>8. Polyolefin fiber filter to reduce contamination.</li> </ol>		

## Product Range

### Individually wrapped in plastic bags

Cat. No.	Volume(mL)	Length(mm)	Sterile	/Pack	/Box	/Case
324003	1	280	YES	1	500	3000
325003	2	280	YES	1	400	2400
326003	5	340	YES	1	200	800
327003	10	339	YES	1	200	800
327303	10 (Short radium carving degree)	244	YES	1	50	500
328003	25	309	YES	1	200	800

**Head Office** No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China  
Tel: +86+ 510-6800 6788      Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
Online: [www.cell-nest.com](http://www.cell-nest.com)

**Oversea** NEST USA (New Jersey/ Phoenix)  
NEST scientific 株式会社 (Yokohama, Japan)  
NEST Scientific Europe B.V (Netherlands)  
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

329003	50	344	YES	1	100	600
--------	----	-----	-----	---	-----	-----

Individually wrapped in paper-plastic bags

324001	1	280	YES	1	500	3000
325001	2	280	YES	1	400	2400
326001	5	340	YES	1	200	800
327001	10	339	YES	1	200	800
328001	25	309	YES	1	200	800
329001	50	344	YES	1	100	600
329501	100	345	YES	1	50	300







**Product Picture**



**Head Office** No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China  
 Tel: +86+ 510-6800 6788 Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
 Online: [www.cell-nest.com](http://www.cell-nest.com)




**Oversea** NEST USA (New Jersey/ Phoenix)  
 NEST scientific 株式会社 (Yokohama, Japan)  
 NEST Scientific Europe B.V (Netherlands)  
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

## Technical Drawing



**Head Office** No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China  
 Tel: +86+ 510-6800 6788 Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
 Online: [www.cell-nest.com](http://www.cell-nest.com)

**Oversea** NEST USA (New Jersey/ Phoenix)  
 NEST scientific 株式会社 (Yokohama, Japan)  
 NEST Scientific Europe B.V (Netherlands)  
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)



**Head Office** No. 530, Xida Road, Meicun Industrial Park, Xinwu District, Wuxi, Jiangsu, China  
 Tel: +86+ 510-6800 6788 Email: [info@nest-wuxi.com](mailto:info@nest-wuxi.com)  
 Online: [www.cell-nest.com](http://www.cell-nest.com)

**Oversea** NEST USA (New Jersey/ Phoenix)  
 NEST scientific 株式会社 (Yokohama, Japan)  
 NEST Scientific Europe B.V (Netherlands)  
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)