

<b>Product Name:</b>	<b>Pasteur Pipette</b>
<b>Product Application:</b>	It is suitable for drawing, transferring and blending a certain volume of liquid. Both sterile and non-sterile Pasteur pipettes are available.
<b>Product Materials:</b>	<b>PP, LDPE</b>
<b>Product Specifications:</b>	<b>Temperature range:</b> stored at room temperature <b>Shelf life:</b> 2 years after date of production
<b>Sterilization:</b>	<b>Yes.</b>
<b>Pyrogens:</b>	Non-Pyrogenic
<b>RNase/DNase:</b>	DNase/RNase free
<b>BSE/TSE:</b>	These products are deemed animal free.
<b>Performance Testing:</b>	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. <b>Appearance inspection:</b> pass <b>Packaging inspection:</b> pass
<b>Features:</b>	<ol style="list-style-type: none"> <li>1. Googlab Pasteur pipettes are made of PP and manufactured in ISO 8 grade cleaning room.</li> <li>2. Good chemical stability.</li> <li>3. Transparent with clear graduations for easy observation.</li> <li>4. The tube wall has good fluidity, good elasticity and easy control. It is suitable for rapid transfer of liquids.</li> <li>5. Pipette can be twisted for use with micro vessels and special-shaped containers.</li> <li>6. High precision and good repeatability.</li> <li>7. Individually packed pasteur pipettes are sterilized by EO.</li> <li>8. DNase/RNase, endotoxin and pyrogenic free.</li> </ol>

## Product Range:

### Pasteur Pipette

SKU	Volume (ml)	Sterile	Scale	Material	Type of Packaging	Pack	Case
GBPWI300	3	No	Yes	PP	Bag	1x500pcs.	500pcs.
GBPII500	3	Yes	Yes	LDPE	Individual	500x1pcs.	500pcs.

Warehouse in Europe

**GenoPlast Biotech S.A.**  
Street. Brzozowa 8  
83-200 Rokocin  
Poland

**E-mail:** office@googlabsci.com


Warehouse in the USA

**GoogLab Scientific USA LLC**  
969 2nd Street SE  
Charlottesville, VA 22902

**E-mail:** office@googlabsci.com

**Technical Drawing:**

**3ml Pasteur pipette**



Warehouse in Europe

**GenoPlast Biotech S.A.**  
Street. Brzozowa 8  
83-200 Rokocin  
Poland

**E-mail:** [office@googlabsci.com](mailto:office@googlabsci.com)

Warehouse in the USA

**GoogLab Scientific USA LLC**  
969 2nd Street SE  
Charlottesville, VA 22902

**E-mail:** [office@googlabsci.com](mailto:office@googlabsci.com)